

GO BIONIC

with the iLet Bionic Pancreas

No more:

- ⊘ Carb counting*
- ⊘ Correction factors
- ⊘ Carb ratios
- ⊘ Pre-set basal rates

or any of the other settings that might be overwhelming about other insulin delivery devices.

The iLet needs only one number - your child's weight.

After startup, the iLet learns about your child - their basal insulin needs, how to correct them to target, and how much insulin they really need for meals — no more guessing!

**For more information
please call Premier Kids Care
at 888-892-9001**



Beta Bionics


betabionics.com

Safety Information - The iLet Bionic Pancreas System is indicated for use by people with type 1 diabetes 6 years of age and older. The iLet Bionic Pancreas requires prescription by a physician. Refer to the iLet Bionic Pancreas System User Guide at www.betabionics.com/resources/user-guides, or for complete safety information including indications, contraindications, warnings, cautions, compatible devices, compatible drugs and instructions, refer to www.betabionics.com/safety. DO NOT start to use the iLet Bionic Pancreas System without adequate training. Incorrect use may result in over-delivery or under-delivery of insulin, which could lead to hypoglycemia or hyperglycemia.

Medical Disclaimer - This handout is for information only and is not a substitute for medical advice and/or services from a healthcare provider. All personal health care decisions and treatment should be discussed with a healthcare provider who is familiar with your individual needs.



*User must be carb aware.
REF: MM000023 Rev A 08/2023



Why Go Bionic?

The **iLet Bionic Pancreas** automatically adjusts to your child's insulin needs so you don't have to.



No carb counting

Instead of carb counting everything on your child's plate, the iLet just needs an estimate of the carbs in the meal - Usual for Me, More or Less.



No corrections

The iLet makes the corrections for you. If your child is trending high, it will automatically give them more insulin to bring them back in range. If they're trending low, it will reduce the current insulin dosage.



iLet algorithms

The "Bionic" part of the iLet that communicates with your child's CGM and makes all of the insulin dosing decisions—no carb counting, no correction factors, no adjustments.



No pre-set basal rate

The iLet learns your child's needs and automatically adjusts their basal rate on the go.



Simple setup

No settings are required to be entered or adjusted—the iLet system only needs your child's weight to get started.



Talks to the Dexcom CGM

The iLet communicates directly with your child's Dexcom G6 to know where their glucose is headed.



Small but mighty

Small (2.3" x 3.5" x 0.6"), lightweight and easy to put in a pocket or clip onto clothing.



Water resistant

If the iLet is dropped into water—no problem! It's protected in up to 12 feet of water for up to 30 minutes.